

ANAEROBIC ADHESIVE

EVOBOND ANAEROBIC ADHESIVE

High Performance anaerobic adhesive is a single component adhesive and cures in the absence of air. EVOBOND anaerobic also cures to form a rigid and durable film to prevent loosening out of vibration and leakage out of corrosion.

Main Features of EVOBOND:

- * THREAD LOCKING SERIES
- * SEALING SERIES
- * RETAINING SERIES




THREAD LOCKING SERIES

Grade	Color	Viscosity(CPS) GB/T2794-1995	Density (g/cm ³)	Max. Gap Fill (mm)	Breakaway Torque (N.m)	Breakloose Torque (N.m)	Work Temp. (°C)	Fixing time (min)	Full Curing Time (h)
AA1222	Purple	1200/thixotropic	1.05	0.12	5.0	3.2	-60 ~ 150	10 ~ 15	24
AA1243	blue	2000/thixotropic	1.09	0.15	14.0	8.0	-60 ~ 150	10 ~ 20	24
AA1270	green	600	1.10	0.13	30.0	25.0	-60 ~ 150	10 ~ 20	24

FLANGE SEALING SERIES


Grade	Color	Viscosity(CPS) GB/T2794-1995	Density (g/cm ³)	Gap Fill (mm)	Work Temp. (°C)	Fixing time (min)	Full Curing Time (h)
AA1510	Pink	40000~140000	1.10	0.25	-60~204	5~10	24
AA1515	Mauve	150000~375000	1.10	0.25	-60~150	5~10	18
AA1518	Red	500000~1000000	1.10	0.25	-60~150	5~10	24



EVOBOND UV CURING ADHESIVE

Professional adhesives for bonding all kinds of transparent substance such as Glass, Crystal, Metal, and so on.

Grade	CPS	Shore	Substrates
T-8086	7600	D71	Glass to Glass
T-8084	20000	D71	Glass to Glass
T-8390	1000	D73	Glass to Metal
T-8085	6300	D72	Glass to Metal
T-8103	1200	D52	Acrylic to Acrylic
T-8323	4000	D68	Acrylic to Acrylic
T-8246	100	D50	PVC to Acrylic
T-8333	9000	D75	Glass to Plastic



CYANOACRYLATE ADHESIVE

CONSUMER PRODUCT



SG-2001F SG-02X EB-541 SG-0301 SG-0301+AD-1 SG-0301 GEL SG-0301 GEL+ NG-10 NG-BO series (5g, 15g)



EV-5 N.T.20g ES-805 N.T.24g (SYRINGE PACKING) ES-8048 N.T.8g (SYRINGE PACKING) PT-2S SG-0112 SG-0220 SG-0312 SG-0312 GEL

INDUSTRY PRODUCT



SH20-20gA LH50-50gA SB20-20gF SB20-20gA



SH20-20gCF SH20-20gF SA50 YB500

TONG SHEN ENTERPRISE CO., LTD.

TAIWAN:
NO.116, LANE 130th, DING-TSO RD., LIN-YUAN DISTRICT,
KAOHSIUNG CITY, TAIWAN
TEL: 886-7-642-0701~5 FAX: 886-7-642-0700
E-mail: evotite@tongshen.com.tw
Http://www.tongshen.com.tw

